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
Laboratory Report

EliA IgA

Run number: 120416-590
Instrument name: 250_1
Instrument type: Phadia 250
Status: Ready
Started: 16-04-2012 10:41:56
Ended: 16-04-2012 14:54:57
Approved: Not approved
Calibrators in run: No

Calibration Information

Measured time: 22-03-2012 13:34:59
Conjugate lot: BVRCD
Acceptance: OK
Active: Yes
Calculation type: Rodbard 4-par
Using zero-point: 1,4

Calibrator	Conc	Measured	Resp	Calc	Judgement	Auto Resp	Edited	Remarks
CAL-0	0,00	Mean	32	0,00				
		13:44:58	32	0,00	OK			
		13:45:58	32	0,00	OK			
		%CV	0,8	0,0				
CAL-0.3	0,30	Mean	190	0,30				
		13:42:58	189	0,30	OK			
		13:43:58	192	0,30	OK			
		%CV	1,2	1,4				
CAL-1.5	1,50	Mean	869	1,51				
		13:40:58	865	1,50	OK			
		13:41:58	872	1,52	OK			
		%CV	0,6	0,6				
CAL-5	5,00	Mean	2630	4,94				
		13:38:59	2627	4,93	OK			
		13:39:59	2633	4,95	OK			
		%CV	0,2	0,2				
CAL-15	15,0	Mean	6883	15,1				
		13:36:59	6803	14,9	OK			
		13:37:59	6963	15,4	OK			
		%CV	1,6	2,1				
CAL-80	80,0	Mean	19467	79,8				
		13:34:59	19403	79,2	OK			
		13:35:59	19531	80,4	OK			
		%CV	0,5	1,1				

Curve data	Value	Logged mean	Difference
ED-20	7,66	6,98	+9,84
ED-50	24,0	21,7	+10,75
ED-80	50,5	46,5	+8,77
Slope	288	292	-1,37

Number of curves included in mean: 106

Curve Controls

Identity	Expected	Measured	Resp	Conc	Judgement	Remarks
CC-1	5,00	Mean	2897	5,5		
		10:41:56	2841	5,4		
		10:42:56	2954	5,6	OK	
		%CV	2,8	3,1	OK	

Quality Controls

Identity	Test	Dil	Measured	Resp	Conc/Quot	Lot	Judgement	Remarks
Celi Pos	Acy	100	10:44:56	14131	77 U/ml	C9XB5	OK	
EliA Neg	Acy	100	10:43:56	54	0,1 U/ml	C9ZB2	OK	


Laboratory Report

EliA IgA

Run number: 120418-594
Instrument name: 250_1
Instrument type: Phadia 250
Status: Ready
Started: 18-04-2012 09:59:37
Ended: 18-04-2012 13:30:40
Approved: Approved
Calibrators in run: No

Calibration Information

Measured time: 22-03-2012 13:34:59
Conjugate lot: BVRCD
Acceptance: OK
Active: Yes
Calculation type: Rodbard 4-par
Using zero-point: 1,4

Calibrator	Conc	Measured	Resp	Calc	Judgement	Auto Resp	Edited	Remarks
CAL-0	0,00	Mean	32	0,00				
		13:44:58	32	0,00	OK			
		13:45:58	32	0,00	OK			
		%CV	0,8	0,0				
CAL-0.3	0,30	Mean	190	0,30				
		13:42:58	189	0,30	OK			
		13:43:58	192	0,30	OK			
		%CV	1,2	1,4				
CAL-1.5	1,50	Mean	869	1,51				
		13:40:58	865	1,50	OK			
		13:41:58	872	1,52	OK			
		%CV	0,6	0,6				
CAL-5	5,00	Mean	2630	4,94				
		13:38:59	2627	4,93	OK			
		13:39:59	2633	4,95	OK			
		%CV	0,2	0,2				
CAL-15	15,0	Mean	6883	15,1				
		13:36:59	6803	14,9	OK			
		13:37:59	6963	15,4	OK			
		%CV	1,6	2,1				
CAL-80	80,0	Mean	19467	79,8				
		13:34:59	19403	79,2	OK			
		13:35:59	19531	80,4	OK			
		%CV	0,5	1,1				

Curve data	Value	Logged mean	Difference
ED-20	7,66	6,98	+9,84
ED-50	24,0	21,7	+10,75
ED-80	50,5	46,5	+8,77
Slope	288	292	-1,37

Number of curves included in mean: 106

Curve Controls

Identity	Expected	Measured	Resp	Conc	Judgement	Remarks
CC-1	5,00	Mean	2663	5,0		
		09:59:37	2694	5,1	OK	
		10:00:37	2631	4,9	OK	
		%CV	1,7	1,8		

Quality Controls

Identity	Test	Dil	Measured	Resp	Conc/Quot	Lot	Judgement	Remarks
Celi Pos	Acy	100	10:02:37	13964	75 U/ml	C9XB5	OK	
EliA Neg	Acy	100	10:01:37	52	0,1 U/ml	C9ZB2	OK	

Sample Report

Sample ID 0001

Sample date: 16-04-2012
Sample status: Reported
Pre-dilution factor: 1
Rack ID: 0050
Source: Local

Sample type: Patient
Tube type: Normal
Priority: Low
Rack pos: 1

Patient name/address:

Requestor name/address:

Patient ID:
Birth date:
Sex:
Diagnosis:

Request ID: 0001
Request date: 16-04-2012
Requestor ID:

Requested test results

EliA IgA

<u>Test name</u>	<u>Test long name</u>	<u>Conc</u>	<u>Class</u>	<u>Quotient</u>	<u>Cut-off</u>	<u>Cut-off 2</u>
tga	EliA Celikey IgA	0,1 U/ml	Negative			

Sample Report

Sample ID 0002

Sample date: 16-04-2012
Sample status: Reported
Pre-dilution factor: 1
Rack ID: 0050
Source: Local

Sample type: Patient
Tube type: Normal
Priority: Low
Rack pos: 2

Patient name/address:

Requestor name/address:

Patient ID:
Birth date:
Sex:
Diagnosis:

Request ID: 0002
Request date: 16-04-2012
Requestor ID:

Requested test results

EliA IgA

<u>Test name</u>	<u>Test long name</u>	<u>Conc</u>	<u>Class</u>	<u>Quotient</u>	<u>Cut-off</u>	<u>Cut-off 2</u>
tga	EliA Celikey IgA	0,5 U/ml	Negative			

Sample Report

Sample ID 0003

Sample date: 16-04-2012
Sample status: Reported
Pre-dilution factor: 1
Rack ID: 0050
Source: Local

Sample type: Patient
Tube type: Normal
Priority: Low
Rack pos: 3

Patient name/address:

Requestor name/address:

Patient ID:
Birth date:
Sex:
Diagnosis:

Request ID: 0003
Request date: 16-04-2012
Requestor ID:

Requested test results

EliA IgA

<u>Test name</u>	<u>Test long name</u>	<u>Conc</u>	<u>Class</u>	<u>Quotient</u>	<u>Cut-off</u>	<u>Cut-off 2</u>
tga	EliA Celikey IgA	0,7 U/ml	Negative			

Sample Report

Sample ID 0004

Sample date: 16-04-2012
Sample status: Reported
Pre-dilution factor: 1
Rack ID: 0050
Source: Local

Sample type: Patient
Tube type: Normal
Priority: Low
Rack pos: 4

Patient name/address:

Requestor name/address:

Patient ID:
Birth date:
Sex:
Diagnosis:

Request ID: 0004
Request date: 16-04-2012
Requestor ID:

Requested test results

EliA IgA

<u>Test name</u>	<u>Test long name</u>	<u>Conc</u>	<u>Class</u>	<u>Quotient</u>	<u>Cut-off</u>	<u>Cut-off 2</u>
tga	EliA Celikey IgA	0,4 U/ml	Negative			

Sample Report

Sample ID 0005

Sample date: 16-04-2012
Sample status: Reported
Pre-dilution factor: 1
Rack ID: 0050
Source: Local

Sample type: Patient
Tube type: Normal
Priority: Low
Rack pos: 5

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Ver. 05 2012.03.01

TABLE OF CONTENTS

1. WARNINGS AND PRECAUTIONS	11
2. DEFINITIONS	11
3. IDENTIFICATION	12
4. PRODUCT SPECIFICATIONS	12
5. TRANSPORTATION AND STORAGE	13
6. PRINCIPAL	13
7. INSTALLING YOUR DEVICE	13
8. PROCEDURE	13
9. CLEANING	17
10. MAINTENANCE	18
11. OUT OF SERVICE	18
12. CE DECLARATION OF CONFORMITY	19

This document is conform to the original version, in Italian and/or English, which is the standard reference.

1. WARNINGS AND PRECAUTIONS

Before using the instrument, carefully read the present manual and the instructions for use of the reagent kit used. Follow all indications on the labels and the instructions provided by the manufacturer. If you have any doubts in regards to the procedure, kindly contact BioDiagene's technical assistance or your local distributor. BioDiagene recommends that all users comply with national standards for health and safety, for example precautions for protection. These protective precautions may include, goggles, gloves and laboratory clothing worn during use or maintenance of this or other laboratory instruments.

If the instrument is used for purposes other than those specified by the manufacturer, safety may be jeopardized.

BioDiagene is committed to meet the specific requirements regarding an "ecological" disposal of its products. The goal is to reduce waste incurred by electrical and electronic equipment.

BioDiagene is aware of the benefits that arise from the possible reuse, treatment, recycle or recovery of devices to minimize the harmful substances emitted into the environment.

Users are responsible for guaranteeing that the devices that contain the symbol indicated on the side reference are not disposed in local waste bins, but that local authority be contacted for correct disposal.

The information contained in this document is subject to change without warning.

There are no user serviceable parts inside this unit. Do not attempt any repairs yourself. If you require technical assistance, please contact BioDiagene or your local distributor.

When you receive the device and during unpacking, carefully check your device for any damage that any have occurred during shipping.

If damage occurs, do not plug the unit in, please contact BioDiagene or your authorized dealer as soon as possible.

2. DEFINIZIONI

Caution, read carefully.

Caution, risk of electric shock.

The product should not be disposed in your local waste bin.

For other warning symbols see the UNI CEI EN 980:2009 directive and subsequent updates.

3. IDENTIFICATION

BioRun is a diagnostic in vitro device (IVD) manufactured by BioDiagene S.r.l., containing an internal firmware (not user-upgradable).

The present manual refers to the following models/versions:

- model BDS 200;
- firmware version 1.xx.xxx.006.

4. PRODUCT SPECIFICATIONS

TECHNICAL CHARACTERISTICS

Power supply: 100 ÷ 240VAC, 50/60/400Hz.

Nominal power: 15 VA.

Fuse: n° 1 – 1,5A – 5x20 mm.

The internal fuse of the device can be substituted with the same fuse type having the same dimensions.

The fuse must be substituted by qualified personnel or by contacting the manufacture's service assistance.

Weight: 4 kg.

Dimensions: 1811 x 275w x 170h mm.

Appliance class*: II.

* in accordance with the Low Voltage Directive 2006/95/EEC.

Link to power supply: connector IEC-320 C14.

External connection: the device has a USB connector (type B).

INFORMATION FOR SECURITY PURPOSES

Control that the electrical system respects the corresponding directive for electrical installations.

ENVIRONMENTAL CONDITIONS

Environmental temperature for correct functioning: T_{max}= 35°C; T_{min}= 0°C.

GENERAL WARNINGS

Keep away from heat sources.

For any technical support please refer to INSTALLING YOUR DEVICE and MAINTENANCE section of the present manual.

SYSTEM FAILURE/DAMAGE

If system failure or damage occurs, contact the manufacture or authorized local distributor: BioDiagene S.r.L., via Aquileia, 34, 90144 Palermo, Italy Tel. +39 0916932138 e-mail info@biodiagene.com.

5. TRANSPORTATION AND STORAGE

There are no specific instructions for the transportation and storage of the device. The use of the original packaging is recommended during transportation.

6. PRINCIPAL

BioRun is a fluorescent reader, able to solicit and measure the emission of fluorophore probes. **The instrument works exclusively with the diagnostic kits produced by BioDiagene.**

The following lists the wavelengths read by the device (excitation/emission):

- 525nm/588nm (HEX);
- 460nm/520nm (FAM).

The biological samples that will be tested by the device are placed in the device using the PCR tube rack containing 12 wells for 0.2 ml PCR tubes (1 strip).

The tests results are visualized on the touchscreen screen display of the device and stored (up to a maximum of 256 analyses) in its non volatile internal memory. The results may also be transferred and stored to a personal computer using the interface software and a USB.

7. INSTALLING YOUR DEVICE

CONNECTION WITH THE EXTERNAL ELECTRICAL SYSTEM

Before starting with the device installation, assure that the characteristics of the electrical system used for power supply are compatible with the technical characteristics of the device used (see section 4 Product Specifications).

Use the power cord supplied with the device; do not use the power cord if it shows signs of lacerations or abrasions on the surface or ends.

8. PROCEDURE

START TURN-ON

On the back of the device you will find a power switch, located near the power cord connector. It is recommended that you position the switch in the off position, switching it on only when necessary.

When the instrument is switched on, the interactions between the operator and the device occur as follows:

- icons found on the touchscreen display;
- RFID badge supplied with the BioDiagene diagnostic kits.

When the instrument is turned on, the following image is displayed on the screen for about 5 seconds.

Image 1 – Start display

Touch the screen and the main menu will appear. Here you can access the functions described in the following paragraphs.

Image 2 – Main menu

In each screen the icons, HOME and BACK are always displayed, respectively in the bottom left and right hand side of the screen. The Home icon is represented by a rectangle containing BioDiagene's logo. The BACK icon is represented by an arrow pointing to the left hand side of the screen. Touching the HOME button takes you back to the main menu, whereas the BACK icon takes you back to the previous screen.

Image 3 – HOME icon (left) and BACK icon (right)

PROCEDURE FOR RUNNING ANALYSES

An analysis is activated by touching the icon ANALYSIS of the device, located on the left hand side of the main menu (*Image 2*). Once the procedure is activated, the device must identify the production lot of the diagnostic kit being used for the analysis. In order to identify the lot simply place the badge, found in the diagnostic kit being used, in the area under the touchscreen, as illustrated on the display of the device (*Image 4*).

The device, after 3 seconds, recognizes the kit and production lot used. This procedure is verified with the screen represented in *Image 5*.

Next, touch the OK icon to add patient information (this step is optional).

Based on the diagnostic kit used, one analysis session is equal to a maximum of 12 biological samples tested. The device also guides you in adding patient information, necessary for the specific diagnostic kit used, with the help of images that appear on the display. This allows for correct data entry that corresponds to the position of the tubes placed in the device.

Nevertheless, for each patient the device automatically associates a specific ID, visualized on the display.

For more information please see the instructions for use (IFU) of the diagnostic kit being used.

After you have completed patient data entry (optional), you can place the PCR tubes into the device. By touching the icon OPEN (located on the bottom left hand side of the screen), the PCR tube rack is extracted from the device. Next, place the PCR tubes in the rack, touch the icon START and the analysis is activated. After about 30 seconds the reading and interpretation of the results is completed and the results are automatically visualized on the display. Finally, remove the samples from within the device by touching the icon OPEN. The PCR tube rack is extracted and you can easily remove the samples from the rack.

Note: for the handling and disposal of the samples, please see the instructions for use of the BioDiagene diagnostic kit used.

Image 4 – Identification of the diagnostic kit


Image 5 – Badge recognition

Image 6 – OPEN and START icon
(located in the bottom right hand side of the display)

ARCHIVE

By touching the icon ARCHIVE, found on the main menu of the device, you can access the list of the last 256 analyses run. For each test, you can visualize the ID automatically associated to the sample and the patient data for each sample.

Using the icons in the form of an arrow, one facing upward and the other facing down, located on the bottom right hand side of the screen, you can scroll through the last 256 results stored in the device.

By touching ID or surname on the screen you can visualize the results of the patient selected. For more information in regards to the visualization of the results please see the instructions for use (IFU) of the diagnostic kit being used.


Image 7 – display of the last stored analyses

DATA TRANSFER TO PC

Using the software iBio, users are able to view and save all the data memorized onto BioRun to their PC (See the iBio IFU).

OPTIONS MENU

To enter the Option menu (Image 8 – Option menu) touch the icon OPTIONS located on the main menu (see Image 2 – Main menu).

By entering the Option menu you can access the following functions:

- 1.Process Control;
- 2.Contamination Control;
- 3.Checking Instrument;
- 4.Operator.

Image 8 – Option menu

The functions Process Control, Contamination Control and Checking Instrument follow the same procedure:

- touch the function desired;
- touch the icon OPEN (located on the bottom right hand side of the screen);
- allow the PCR tube rack to be extracted from within the device;
- insert the tube containing the required reagents (for detailed information please see the instructions for use (IFU) of the diagnostic kit being used);
- touch the START icon (bottom right hand side);
- place the badge in the area under the touchscreen;
- touch the OK icon;
- wait for the results of the procedure;
- touch the OPEN icon (bottom right hand side);
- remove the tube from the PCR tube rack.

Touching the screen in correspondence to the function Operator brings you to screen indicated in *Image 9*, where you can enter the user's name.

When you have finished entering the requested information, touch the OK icon (bottom right hand side) to confirm. By touching the CANCEL icon (bottom left hand side) the operation is cancelled and the display returns to the Options menu.

Image 9 –Adding operator information

PARTICULAR CASE FINDINGS

If an error occurs while the instrument is working verify table "BIORUN ERROR CODES" (page 18). If the error persists contact the manufacturer or the distributor with the error code.

If the power supply of the device is interrupted or if there is a failure of the analyses, it is recommended to restart the entire analyses procedure. The data pertaining to the last 256 tests stored in the system will not be lost.

9. CLEANING

When cleaning the external areas of the device:

- do not use corrosive products that may damage the plastic components of the device and/or power cords at the expense of electrical safety;
- after cleaning, the surface of the device should be thoroughly dried with a soft cloth;
- do not use detergents that contain alcohol. Dry the polished surface of the instrument thoroughly using a soft cloth in order to avoid wear-down of the polish.

You **MUST** clean the PCR tube rack periodically (with intervals that do not exceed 1 week).

It is strongly recommended that cleaning be performed before running any analyses because contamination may rarely occur in the PCR tube rack due to incorrect handling by the user (PCR strips accidentally open while inside the instrument, placing the PCR strips into the device with the caps slightly raised, etc.).

You can clean the PCR tube rack using a solution composed of 1% bleach and apply to the device using cotton swabs. The procedure should be performed when the PCR tube rack is extracted from the device, that is before placing the tubes in the rack (aka before running the analyses) or after having removed the tubes from the rack, at the conclusion of the test.

The use of other substances other than those specified in the instruction manual for internal cleaning are not permitted.

The use of liquid substances (example spray dispensers) inside the device are not permitted.

10. MAINTENANCE

Checking Instrument

It is necessary to periodically verify that the instrument is functioning correctly by using the kit Checking Instrument (cod. BDF304).

To perform this activity please refer to section 8 Procedure, paragraph Options Menu of the present user manual and the instructions for use of the Checking Instrument kit.

It is recommended to check the instrument every 6 months.

There are no user serviceable parts inside this unit. Maintenance and repair services are authorized only by BioDiagene personnel and authorized dealers. Do not attempt any repairs yourself. In the unlikely event your device may require service, please contact BioDiagene or other authorized dealer.

The device should not be opened by the user. Only the service assistance of the manufacture is permitted to open the device.

11. OUT OF SERVICE

The device is disposed of according to national legislation. No preliminary operations are necessary for disposal.

Do not dispose the device in your local bin.

BIORUN ERROR CODES

Error 7	Unplug the USB cord if plugged to the computer. Restart BioRun by using the power switch. If the error code persists contact the manufacturer
Error 8	Restart BioRun by using the power switch. If the error code persists contact the manufacturer
Error 9	Verify that there are no obstacles that prevent the movement of the PCR tube tray. Restart BioRun by using the power switch. If the error code persists contact the manufacturer

BioDiogene SRL
Via Aquileia, 34/A
90144 Palermo

DECLARES UNDER SOLE RESPONSIBILITY
THAT THE DEVICE FOR IN VITRO DIAGNOSTICS

BioRun (BDS 200)

IS COMPLAINT AS INDICATED IN:

LOW VOLTAGE DIRECTIVE 2006/95/EEC

(Applicable norms: EN 61010-1 / EN 61010-2-081 / EN 61010-2-101)

ELECTROMAGNETIC COMPATIBILITY 89/336/CEE E (AND UPDATES)

(Applicable norms: EN 61326-1 / EN 61326-2-6)

DECLARES THAT THE SYSTEM IS COMPLAINT WITH THE NORMS OF THE CE MARK

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Ver. 02 2012.03.01

INDEX

Application	10
Test Principle	10
Kit Components and Preparation of Reagents	10
Material and Instruments Not Supplied	11
Performance Criteria	11
Clinical Sample	11
Before Starting	12
Warnings and Precautions	12
Procedure	13
Interpretation of Results	16
Validation	16
Troubleshooting	17
References	42

APPLICATION

Kit for the identification of samples susceptible to Celiac Disease performed in one PCR reaction.

TEST PRINCIPLE

Celiac Gene Screen foresees a preliminary lysis of blood sample (EDTA-tube), DNA amplification (1 reaction/test) and fluorescence detection using BioRun Reader (BioDiagene Cat. No. BDS 200).

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EXTRACTION BUFFER	Cat. No. ZBECGS	Blue cap
A bottle containing a ready-to-use extraction buffer. Usage: after collecting the required amount of the solution, tighten and secure the cap, storing at the indicated temperature.		
Storage: +4°C	Stability: 1 year	Stability after opening: 1 year

PRIMER MIX	Cat. No. ZPCGS	Natural 12 tube PCR strip
0,2 ml PCR tubes containing dried primers/probes for a housekeeping gene (internal control) and the target alleles. Usage: open the pouch containing the PCR tubes, collect all necessary tubes and seal with the caps. Replace additional tubes back into their pouch sealed with their corresponding caps and store at the indicated temperature. Attention: do not expose the strips to light for long periods of time.		
Storage: +4°C	Stability: 1 year	Stability after opening: 1 year

TAQ MIX	Cat. No. ZTAQCS	Green cap
2 tubes containing a ready-to-use Mix consisting of Taq Polymerase, MgCl ₂ , dNTPs and buffer. Usage: ready-to-use Mix. Immediately after use, seal with the caps and store the tubes at the indicated temperature.		
Storage: +4°C	Stability: 1 year	Stability after opening: 1 year

PRIMER MIX FOR CONTROL CELIAC GENE	Cat. No. ZCTRCG	Natural 6 tube PCR strip
0,2 ml PCR tubes containing a dried primers/probes mix for a housekeeping gene. Use 1 tube per sample as indicated in "Procedure". Usage: open the pouch containing the PCR tubes, collect all necessary tubes and seal with the caps. Replace additional tubes back into their pouch sealed with their corresponding caps and store at the indicated temperature. Attention: do not expose the strips to light for long periods of time.		
Storage: +4°C	Stability: 1 year	Stability after opening: 1 year

WAX	Cat. No. ZCECGS	White cap
Bottle containing a ready-to-use wax for PCR. Usage: before using allow the wax to reach room temperature. Immediately after use store the bottle at the indicated temperature.		
Storage: +4°C/RT	Stability: 1 year	Stability after opening: 1 year

CAPS FOR 0.2 ML TUBES	Cat. No. ZTAPCGS	Natural caps
Flat natural caps for 0.2 ml tubes in strips of 12 and 6 caps. Usage: use to seal the tubes containing the Primer Mix in strips of 12 and 6 tubes.		
Storage: +4°C/RT	Stability: NA	Stability after opening: NA

BADGE CELIAC GENE SCREEN	Cat. No. ZBADC GS	Natural caps
Badge containing information in regards to the lot number and type of kit used. Usage: approach the badge to the instrument in order for BioRun to read and recognize the lot and application used.		
Storage: +4°C/RT	Stability: NA	Stability after opening: NA

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- Thermal Cycler for 0,2 ml tubes with heated lid
- 10 µl and 200 µl adjustable pipettes
- 10 µl and 200 µl pipette tips with filter
- Disposable gloves
- Sterile 0,5 or 1,5 ml tubes
- Vortex
- Minicentrifuge for strips of 0,2 ml tubes or a centrifuge with adaptor (6000rpm/2000xg)

PERFORMANCE CRITERIA

Analytical Sensibility: the amplification and identification systems guarantee results if lysing at least 2x10³ cells/ reaction.

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Method restrictions: the kit does not discriminate between the haplotype DQ2 and DQ8.

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- Collect blood sample in EDTA-tubes for blood.

- Use blood sample as indicated in the instructions for the EDTA tubes.
- Store blood EDTA tubes at +4°C.
- Do not use the blood sample when agglutinated.
- Work in sanitized environments.
- Do not use the same consumable materials (tubes, vials) for different samples.
- The temperature in the laboratory should not exceed +25°C.
- Use a calibrated refrigerator to store blood samples.
- Sanitize the refrigerator maintaining it in ideal conditions in order to guarantee the appropriate temperature at all times.
- Adapt the appropriate operative precautions in order to avoid contamination between samples.

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N.B.: The amplification program, when using the Celiac Gene Screen kit is identical to the amplification program for the entire Celiac Gene product line.

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Procedural warnings

- It is necessary to remove the PCR tube tray positioned on the Peltier tray before starting amplification.
- Place the PCR tubes directly on the Peltier thermal cycler tray.
- After use close and seal PCR tubes and reagents, storing reagents at their accurate storage temperature. It is very important not to leave the tubes containing the dried primers opened. Seal caps (ZIAPCGS) onto tubes with dried primers. Caution, a humid environment may inactivate the primers.
- Do not expose the Primer Mix strips to the light for long periods of time.
- Do not use the reagents after the expiration date.
- The lack of maintained and/or inaccurate storage temperature causes failure of the IVD.
- A lack of and/or incorrect maintenance of lab instruments jeopardizes results causing IVD failure.
- Do not change or modify the procedure.
- Do not substitute reagents with other lots or made from other manufactures.
- DO NOT USE THE BADGE CELIAC GENE SCREEN WITH DIFFERENT LOTS.
- Do not dispense and/or transfer solutions from their original bottles into different bottles/vials.
- After you have performed the test do not re-use the PCR tubes: the system is disposable (for single-use only).

- Contact the manufacturer in case the IVD package is damaged.
- Use 2 separate working areas: DNA Area (DNA extraction) and PCR Area (for preparation of the reactions). The DNA area is intended for the extraction of nucleic DNA using dedicated pipettes. The PCR Area is intended exclusively for the preparation of the PCR tubes, preferably under a laminar flow Workstation using aerosol barrier tips. Use devices and lab material at their respective place at all times, never exchanging them. Note: Carry-Over contamination due to prior amplification greatly jeopardizes results.
- Avoid bacterial contamination and/or cross-contamination of reagents.
- Do not touch or wet the rim of the vials containing the reagents with tips and pipettes.
- Avoid blowing or speaking over opened vials or tubes.
- The use of regularly calibrated instruments are strongly recommended. Note: use a calibrated thermal cycler only.
- After the PCR program has been completed, do not place the strips, from the thermal cycler, on the same rack that was used for the PCR preparation.
- Two sets of pipettes are necessary for test execution: one set for the DNA extraction phase and the second set when preparing for PCR.
- Do not use direct artificial lighting if you are using a laminar flow Workstation during the preparation of the PCR tubes phase.
- **Read the strips containing the PCR product using BioRun within 1 hour after the PCR program has been completed.**

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- Vortex the tubes containing the blood samples.
- For each blood sample, pipette 200 µl of Extraction buffer into a 1.5 or 0.5 ml sterile tube and add 10 µl of the blood sample.
- Mix or vortex and incubate for 1 minute at room temperature.
- Mix or vortex until reaching a homogeneous color of the lysate.
- **Use 2 µl of lysate for each sample that will be amplified.**

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PCR

- Use one Primer Mix PCR tube (ZPCGS) per sample.
- **Seal the remaining tubes containing the Primer Mix with their corresponding caps (ZTAPCGS) and store at +4°C.**
- Add 18 µl of Taq Mix and 2 µl of DNA to each marked tube. Mix by pipetting. Add 30 µl of wax. Mix by vortexing until the dried primers are completely dissolved. Collect the drops to the bottom of the tubes with a minicentrifuge.
- Store the Taq Mix (tube with green cap) and the wax (tube with white cap) and store at the indicated temperature.
- **Important: Remove the PCR tube tray within the thermal cycler. Do not use at any time!**
- Place the PCR tubes directly onto the Peltier tray in the thermal cycler. Perform the amplification program as previously described.
- Wait until the amplification process is complete (1,5 hours) and, **after the wax reaches room temperature, place the PCR strip in the BioRun instrument in order to read the results. Position the tubes or the strip starting from the first position of the PCR tube rack (internal side). Do not leave empty wells between one tube and another tube.**

HOW TO USE THE BADGE (ZBADCGS) WITH THE INSTRUMENT BIORUN

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HOW TO USE BIORUN TOGETHER WITH THE CELIAC GENE SCREEN KIT

How to start an analysis using BioRun

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How to add patient information

Using Celiac Gene Screen (CGS) you can run a maximum of 12 samples/tube during one work session using BioRun.

Each tube can be associated to a different patient. BioRun uses a graphic interface to help guide the user in adding patient information (optional step) correctly.

When running the Celiac Gene Screen analysis, the display of the device will present an image that represents the PCR tube tray extracted from BioRun with each well position numbered from 1 to 12. On the left hand side of the screen you will see 12 numbered icons with an image of a tube, each associated to a well position of the PCR tube rack.


Image 1 – Celiac Gene Screen tube/patient positioning

Touching the icons takes you to a new display (Image 2), where you can visualize the ID that is automatically assigned by the instrument and the name of the user, if registered.

Image 2 – CGS patient data entry

Image 3 – Icon for the activation of the on-screen keyboard

By touching the Pencil icon (Image 3) you can access the on-screen keyboard and

enter and/or modify patient information, such as name and surname of the patient (Image 4).

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Image 4 – On-screen keyboard

Start reading

After you have added the patient information, from the display represented by Image 1 touch the OPEN icon (bottom right hand side) and the PCR tube tray will be extracted from within the device. **Place the PCR tubes ready to be tested starting from the first position of the PCR tube rack (internal side). Do not leave empty wells between one tube and another tube.**

Touch the START icon (bottom right hand side) to start the analysis. After about 30 seconds the analysis procedure is completed and the results are automatically visualized on the display (Image 5).


Image 5 – Celiac Gene Screen result screen

The results of the analyses of each tube are represented by the following symbols:

+	sample susceptible to Celiac Disease
-	sample not susceptible to Celiac Disease
	PCR tube not present
	Error (repeat analyses)

Note: Errors can occur if the sample was prepared incorrectly or if the operator added patient information for a corresponding well-position where there is no PCR tube placed (tube missing). Finally, after the analysis is completed, you can remove the samples from within the device by touching the icon OPEN (bottom right hand side).

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- Use one sample of DNA extracted in the same work session. Gather one tube of the Primer Mix for Control Celiac Gene (ZCTRCG).
- Add 18 µl of Tax Mix (ZTAQCGS) and 2 µl of DNA. Mix by pipetting. Add 30 µl of wax.
- Mix by vortexing until the dried primers/probes are completely dissolved. Collect the drops to the bottom of the tubes with a minicentrifuge.
- Close the remaining tubes with their corresponding caps. Store at +4°C.
- Perform amplification and activation.
- Touch the icon OPTIONS on the touchscreen display of the BioRun reader.
- Touch the function PROCESS CONTROL.
- Touch the icon OPEN and the PCR tube tray will be extracted.
- Allow the PCR tube rack to be extracted from within the device.
- Insert the Primer Mix for Control Celiac Gene tube on the PCR tube tray.
- Touch the START icon (bottom right hand side).
- Place the badge in the area under the touchscreen.
- Touch the OK icon.
- Wait for the results of the procedure.
- Touch the OPEN icon (bottom right hand side).
- Remove the tube from the PCR tube rack.

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- Gather one tube of the Primer Mix for Control Celiac Gene (ZCTRCG).
- Add 18 µl di Taq Mix (ZTAQCGS). Add 30 µl of wax.
- Mix by vortexing until the dried primers/probes are completely dissolved. Collect the drops to the bottom of the tubes with a minicentrifuge.
- Close the remaining tubes with their corresponding caps (ZTAPCGS). Store at +4°C.
- Perform amplification and activation.
- Touch the icon OPTION on the touchscreen display of the BioRun reader.
- Touch the function CONTAMINATION CONTROL.
- Touch the icon OPEN and the PCR tube tray will be extracted.
- Allow the PCR tube rack to be extracted from within the device.
- Insert the Primer Mix for Control Celiac Gene tube on the PCR tube tray.
- Touch the START icon (bottom right hand side).
- Place the badge in the area under the touchscreen.
- Touch the OK icon.
- Wait for the results of the procedure.
- Touch the OPEN icon (bottom right hand side).
- Remove the tube from the PCR tube rack.

INTERPRETATION OF RESULTS

The interpretation of the results is run automatically by the BioRun reader (BioDiagene Cat. No. BDS 200) as follows:

POSITIVE (symbol "+") : sample susceptible to Celiac Disease

NEGATIVE (symbol "-"): sample not susceptible to Celiac Disease

REPEAT (symbol ): The test must be repeated starting from the DNA extraction from the biological sample.

If the word "CONTAMINATION" appears when running the contamination control, this means that there is contamination with genomic DNA. If this situation occurs please read the procedural warnings carefully.

The samples that result susceptible to Celiac Disease can be typed using the Celiac Gene Typing kit (Cat. No. BDF 302) and Celiac Alleles (Cat. No. BDF 303) for the identification of the alleles associated with Celiac Disease.

VALIDATION

The test is valid when the symbols "+" or "-" appear on the BioRun touchscreen.

If the symbol  appears, the test is not valid and the sample must be retested, starting from DNA extraction.

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TROUBLESHOOTING

PROBLEMA	CAUSA	SOLUZIONE
 The symbol appears on the touchscreen	DNA absent	Repeat DNA extraction. If the problem continues repeat blood withdrawal
	Blood collected in the wrong tube (i.e. collected in a tube containing heparin)	Repeat blood sample withdrawal and store the sample in an EDTA-tube
	Agglutinated or coagulated blood sample	Do not use; repeat blood sample withdrawal
	Instruments (thermal cycler) not calibrated	Use and/or calibrate the lab instruments used for performing the tests
	PCR session failed due to improper use of the thermal cycler	Place PCR tubes directly onto thermal cycler plate, do not use PCR tube rack. Control amplification program
	Evaporation in reaction tubes, due to water loss and change in the concentration of the PCR reagents	Tighten and secure tubes with caps - vortex and centrifuge to collect drops before using thermal cycler
"ERROR 2" appears	The fluorescence values for the identification of the internal control and for the identification of the target allele are low. The primer mix located on the bottom of the tube has not been sufficiently dissolved	Repeat the test and vortex each tube of the strip adequately
"ERROR 4" appears	The fluorescence value for the identification of the target allele is in the uncertain range. The primer mix located on the bottom of the tube has not been sufficiently dissolved	Repeat the test and vortex each tube of the strip adequately
"ERROR 5" appears	No tube detected in the well-position	Verify that the data entry corresponds correctly to the position of the tubes placed in the device.
The word "NEGATIVE" appears in the process control	PCR session failed	You must repeat the test, starting from DNA lysis
The word "CONTAMINATION" appears in the contamination control	Contamination due to amplified products from previous work sessions	Sanitize the work area and instruments

Note: if an error message is displayed on the screen of the BioRun reader kindly contact the manufacture and communicate the error code.

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Ver. 02 2012.03.01

INDEX

Application	10
Test Principle	10
Kit Components and Preparation of Reagents	10
Material and Instruments Not Supplied	11
Performance Criteria	12
Clinical Sample	12
Before Starting	12
Warnings and Precautions	12
Procedure	13
Interpretation of Results	16
Validation	17
Troubleshooting	17
References	42

APPLICATION

Kit for the detection of HLA class II alleles associated with Celiac Disease and for the direct identification of samples homozygous for DQB1*02 allele. The alleles detected are: HLA II DQA1*0201, DQA1*03, DQA1*05, DQB1*02, DQB1*0301/0304, DQB1*0302/0305. Celiac Gene Alleles identifies the risk of developing Celiac Disease based on the targeted alleles.

TEST PRINCIPLE

Celiac Gene Alleles foresees a preliminary lysis of blood sample, DNA amplification (PCR) and fluorescence detection using BioRun Reader. 7 PCR reactions (1 strip) are necessary to carry out one test.

KIT COMPONENTS AND PREPARATION OF REAGENTS

The kit is composed of the following reagents:

EXTRACTION BUFFER	Cat. No. ZBECGA	Blue cap
A bottle containing a ready-to-use extraction buffer. Usage: after collecting the required amount of the solution, tighten and secure the cap, storing at the indicated temperature.		
Storage: +4°C	Stability: 1 year	Stability after opening: 1 year

PRIMER MIX	Cat. No. ZPCGA	Natural 7 tube PCR strip
0,2 ml PCR tubes containing dried primers/probes for a housekeeping gene (internal control) and the target alleles. Usage: open the pouch containing the PCR tubes, collect all necessary tubes and seal with the caps. Replace additional tubes back into their pouch and store at the indicated temperature. Attention: do not expose the strips to light for long periods of time.		
Storage: +4°C	Stability: 1 year	Stability after opening: 1 year

TAQ MIX	Cat. No. ZTAQCGA	
20 tubes containing a ready-to-use Mix consisting of Taq Polymerase, dNTPs and buffer. Usage: ready-to-use Mix. Immediately after use, seal with the caps and store the tubes at the indicated temperature.		
Storage: +4°C	Stability: 1 year	Stability after opening: 1 year

PRIMER MIX FOR CONTROL CELIAC GENE	Cat. No. ZCTRCG	Natural 6 tube PCR strip
0,2 ml PCR tubes containing a dried primers/probes mix for a housekeeping gene. Use 1 tube per sample as indicated in "Procedure". Usage: open the pouch containing the PCR tubes, collect all necessary tubes and seal with the caps. Replace additional tubes back into their pouch and store at the indicated temperature. Attention: do not expose the strips to light for long periods of time.		
Storage: +4°C	Stability: 1 year	Stability after opening: 1 year

TAQ MIX FOR CONTROL CELIAC GENE	Cat. No. ZTAQCCG	Green cap
1 tube containing a ready-to-use Mix consisting of Taq Polymerase, MgCl ₂ , dNTPs and buffer. Usage: ready-to-use Mix. Immediately after use, seal with the cap and store the tube at the indicated temperature.		
Storage: +4°C	Stability: 1 year	Stability after opening: 1 year

WAX	Cat. No. ZCECGA	White cap
Bottle containing a ready-to-use wax for PCR. Usage: before using allow the wax to reach room temperature. Immediately after use tighten and secure the cap and store the bottle at the indicated temperature.		
Storage: +4°C/RT	Stability: 1 year	Stability after opening: 1 year

CAPS FOR 0.2 ML TUBES	Cat. No. ZTAPCGA	Natural caps
Flat natural caps for 0.2 ml tubes in strips of 7 and 6 caps. Usage: use to seal the tubes containing the Primer Mix in strips of 7 and 6 tubes.		
Storage: +4°C/RT	Stability: NA	Stability after opening: NA

BADGE CELIAC GENE ALLELES	Cat. No. ZBADCGA	
Badge containing information in regards to the lot number and type of kit used. Usage: approach the badge to the instrument in order to read and recognize the lot and application used.		
Storage: +4°C/RT	Stability: NA	Stability after opening: NA

REAGENT QUANTITY

Kit	ZBECGA	ZPCGA	ZTAQCGA	ZCTRCG	ZTAQCCG	ZCECGA	ZTAPCGA	ZBADCGA
Cat. No. BDF 303	1 x 7,5 ml	20 strips x 7 tubes	20 x 144 µl	1 strip x 6 tubes	1 x 120 µl	1 x 7,7 ml	20 strips x 7 caps; 1 strip x 6 caps	1 pc

MATERIAL AND INSTRUMENTS NOT SUPPLIED

- Thermal Cycler for 0.2 ml tubes with heated lid
- 10 µl and 200 µl adjustable pipettes
- 10 µl and 200 µl pipette tips with filter
- Disposable gloves
- Sterile 0,5 or 1,5 ml tubes
- Vortex
- Minicentrifuge for strips of 0,2 ml tubes or a centrifuge with adaptor (6000rpm/2000xg)

PERFORMANCE CRITERIA

Analytical Sensibility: the amplification and identification systems guarantee results if lysing at least 2×10^2 cells/reaction.

Analytical Specificity: the composition of the primer/probe mix used in the kit identifies HLA alleles associated with Celiac Disease, based on the latest sequence data and from the analysis of typed DNA samples, extracted from referenced B-lymphoblastoid cell lines, from the IHWG – International Histocompatibility Working Group's Cell and Gene Bank.

Diagnostic Sensibility: a performance study was conducted using 50 DNA samples, typed with a reference IVD (CE marked). The results revealed a diagnostic sensibility of 100%.

Diagnostic Specificity: a performance study was conducted using 50 DNA samples, typed with a reference IVD (CE marked). The results revealed a diagnostic specificity of 100%.

Stability: the kit is stable for 1 year if storage parameters are maintained; stability of an opened kit is equal to 1 year. The kit is stable at room temperature during transportation.

CLINICAL SAMPLE

The kit uses DNA extracted from a blood sample. Blood samples must be collected by well trained personnel, as follows:

- Collect blood sample in EDTA-tubes for blood.
- **Use blood sample as indicated in the instructions for the EDTA tubes.**
- Store blood EDTA tubes at +4°C.
- Do not use the blood sample when agglutinated.
- Work in sanitized environments.
- Do not use the same consumable materials (tubes, vials) for different samples.
- The temperature in the laboratory should not exceed +25°C.
- Use a calibrated refrigerator to store blood samples.
- Sanitize the refrigerator maintaining it in ideal conditions in order to guarantee the appropriate temperature at all times.
- Adapt the appropriate operative precautions in order to avoid contamination between samples.

BEFORE STARTING

1. Allow the tube containing the wax to reach room temperature.
2. Set and memorize the following amplification program to the thermal cycler:

Amplification program

Temperature	Time	Cycles
94°C	2 min	1
94°C	15 sec	40
52°C	60 sec	
15°C	∞	

N.B.: The amplification program, when using the Celiac Gene Alleles kit is identical to the amplification program for the entire Celiac Gene product line.

WARNINGS AND PRECAUTIONS

Storage

- Store the kit at +4°C;
- The kit is stable at room temperature during transportation.

Precautions for personnel

- Wear P.P.E. (Personal Protective Equipment) e.g.: disposable gloves.
- Wash hands thoroughly after handling IVD.
- If reagents come in contact with eyes or skin, immediately rinse abundantly with water.

Procedural warnings

- **Celiac Gene Alleles can not be used for tissue typing.**
- **It is necessary to remove the PCR tube tray positioned on the Peltier tray before starting amplification.**
- Place the PCR tubes directly on the Peltier thermal cycler tray.
- After use close and seal PCR tubes and reagents, storing reagents at their accurate storage temperature. **It is very important not to leave the tubes containing the dried primers opened.** Seal caps (ZIAPCGA) onto tubes with dried primers. Caution, a humid environment may inactivate the primers.
- **Do not expose the Primer Mix strips to the light for long periods of time.**
- Do not use the reagents after the expiration date.
- The lack of maintained and/or inaccurate storage temperature causes failure of the IVD.
- A lack of and/or incorrect maintenance of lab instruments jeopardizes results causing IVD failure.
- Do not change or modify the procedure.
- **Do not substitute reagents with other lots or made from other manufactures.**
- **DO NOT USE THE BADGE CELIAC GENE ALLELES WITH DIFFERENT LOTS.**
- Do not dispense and/or transfer solutions from their original bottles into different bottles/vials.
- Do not expose reagents to strong light.
- After you have performed the test do not re-use the PCR tubes; the system is disposable (for single-use only).
- Contact the manufacturer in case the IVD package is damaged.
- Use 2 separate working areas: DNA Area (DNA extraction) and PCR Area (for preparation of the reactions). The DNA area is intended for the extraction of nucleic DNA using dedicated pipettes. The PCR Area is intended exclusively for the preparation of the PCR tubes, preferably under a laminar flow Workstation using aerosol barrier tips. Use devices and lab material at their respective place at all times, never exchanging them. Note: Carry-Over contamination due to prior amplification greatly jeopardizes results.
- Avoid bacterial contamination and/or cross-contamination of reagents.
- Do not touch or wet the rim of the vials containing the reagents with tips and pipettes.
- Avoid blowing or speaking over opened vials or tubes.
- The use of regularly calibrated instruments are strongly recommended. Note: use a calibrated thermal cycler only.
- After the PCR program has been completed, do not place the strips, from the thermal cycler, on the same rack that was used for the PCR preparation.
- Two sets of pipettes are necessary for test execution: one set for the DNA extraction phase and the second set when preparing PCR tubes.
- Do not use direct artificial lighting if you are using a laminar flow Workstation during the preparation of the PCR tubes phase.
- **Read the strips containing the PCR product using BioRun within 1 hour after the PCR program has been completed.**

PROCEDURE

LYSIS OF BLOOD SAMPLE COLLECTED USING EDTA-TUBES

- Vortex the tubes containing the blood samples.
- For each blood sample, pipette 200 µl of Extraction buffer into a 1.5 or 0.5 ml sterile tube and add 10 µl of the blood sample.
- Mix or vortex and incubate for 1 minute at room temperature.
- Mix or vortex until reaching a homogeneous color of the lysate.
- **Use 16 µl of lysate for each sample that will be amplified.**

Note: All BioDiagene Celiac Gene kits use the same lysis procedure. Therefore one lysate can be used with all the kits of the same product line including the reaction control.

PCR

- Dedicate one strip of 7 tubes of the Primer Mix (ZPCGA) per sample.
- Gather a number of strips equal to the number of samples being tested. Number each strip. The first tube

is marked and contains primers for DQA1*0201 allele.

- Gather the Taq Mix tubes equal to the number of samples being tested. Dedicate one Taq Mix tube per sample.
- For each sample dispense 16 µl of lysate to one Taq Mix tube. Mix by vortexing.
- Pipette 20 µl of this mix solution (Taq Mix and DNA) in each PCR tube for one strip per sample. Add 30 µl of wax. Mix by vortexing until the dried primers/probes are completely dissolved.
- Collect the drops to the bottom of the tubes with a minicentrifuge.
- **Close the remaining tubes and seal the Primer Mix tubes with their corresponding caps (ZTAPCGA). Store at +4°C.**
- **Important: Remove the PCR tube tray within the thermal cycler. Do not use at any time!**
- Place the PCR tubes directly onto the Peltier tray in the thermal cycler. Perform the amplification program as previously described.
- Wait until the amplification process is complete (1,5 hours) and, **after the wax reaches room temperature, place the PCR strip in the BioRun instrument in order to read the results.**
- **Place the PCR tubes, ready to be tested, with the marked tube into the tray starting from the first position of the PCR tube rack (internal side).**

HOW TO USE THE BADGE (ZBADCGA) WITH THE INSTRUMENT BIORUN

The Badge identifies the type of application used and the lot number. The badge must be used before you perform a reading using BioRun.

To activate the analysis procedure it is necessary that the instrument identifies the lot number of the production kit being used. In order to identify the lot simply place the badge, found in the kit being used, in the area under the touchscreen, as illustrated on the display of the device, and wait for about 3 seconds.

Make sure the lot number written on the kit matches the lot number of the badge for every reading/analysis.

HOW TO USE BIORUN TOGETHER WITH THE CELIAC GENE ALLELES KIT

Please see the instructions for use of the BioRun reader.

Read the strips containing the PCR product using BioRun within 1 hour after the PCR program has been completed.

How to add patient information

The Celiac Gene Alleles kit allows you to run 1 analysis per patient sample (7 tubes) using the BioRun reader.

When you start an analysis, a screen appears that allows you to add patient information (including ID number and the name of the user, if registered).


Image 1 – CGA patient data entry

By touching the pencil icon (Image 2) you can access the on-screen keyboard and

Image 2 – Icon for the activation of the on-screen keyboard

enter and/or modify patient information, such as name and surname of the patient (Image 3).

When adding patient information, you can touch the OK icon (bottom right hand side) and CANCEL icon (bottom left hand side) to confirm or cancel the operation.

Start reading

After you have added the patient information, from the display represented by Image 1 touch the OPEN icon (bottom right hand side) and the PCR tube tray will be extracted from within the device. **Place the PCR tubes, ready to be tested, with the marked tube into the tray starting from the first position of the PCR tube rack (internal side).**

Touch the START icon (bottom right hand side) to start the analysis. After about 30 seconds the analysis procedure is completed and the results are automatically visualized on the display (Image 4).

The Celiac Gene Alleles results are visualized in two different screens; the first screen lists all the identified alleles, the second screen is accessed by touching the RESULT icon and lists the serologic haplotype given the identified alleles.

Finally, after the analysis is completed remove the PCR tubes from within the device, by touching the OPEN icon (bottom right hand side).

Image 3 – On-screen keyboard

Image 4 – CGA Result interpretation

PRIMER MIX FOR CONTROL CELIAC GENE (COD. ZCTRGC)

Follow the instructions below when running the **process control**:

- Use one sample of DNA extracted in the same work session. Gather one tube of the Primer Mix for Control Celiac Gene (ZCTRGC).
- Add 18 µl of Tax Mix for Control Celiac Gene (ZTAQCCG) and 2 µl of DNA. Mix by pipetting. Add 30 µl of wax.
- Mix by vortexing until the dried primers/probes are completely dissolved. Collect the drops to the bottom of the tubes with a minicentrifuge.
- Close the remaining tubes with their corresponding caps. Store at +4°C.
- Perform amplification.
- Touch the icon OPTIONS on the touchscreen display of the BioRun reader.
- Touch the function PROCESS CONTROL.
- Touch the icon OPEN and the PCR tube tray will be extracted.
- Allow the PCR tube rack to be extracted from within the device.
- Insert the Primer Mix for Control Celiac Gene tube on the PCR tube tray.
- Touch the START icon (bottom right hand side).
- Place the badge in the area under the touchscreen.
- Touch the OK icon.
- Wait for the results of the procedure.
- Touch the OPEN icon (bottom right hand side).
- Remove the tube from the PCR tube rack.

Follow the instructions below when running the **contamination control**:

- Gather one tube of the Primer Mix for Control Celiac Gene (ZCTRGC).
- Add 18 µl di Taq Mix for Control Celiac Gene (ZTAQCCG). Add 30 µl of wax.
- Mix by vortexing until the dried primers/probes are completely dissolved. Collect the drops to the bottom of the tubes with a minicentrifuge.
- Close the remaining tubes with their corresponding caps (ZTAPCA). Store at +4°C.
- Perform amplification.
- Touch the icon OPTION on the touchscreen display of the BioRun reader.
- Touch the function CONTAMINATION CONTROL.
- Touch the icon OPEN and the PCR tube tray will be extracted.
- Allow the PCR tube rack to be extracted from within the device.
- Insert the Primer Mix for Control Celiac Gene tube on the PCR tube tray.
- Touch the START icon (bottom right hand side).
- Place the badge in the area under the touchscreen.
- Touch the OK icon.
- Wait for the results of the procedure.
- Touch the OPEN icon (bottom right hand side).
- Remove the tube from the PCR tube rack.

INTERPRETATION OF RESULTS

The interpretation of the results is run automatically by the BioRun reader (BioDiagene Cat. No. BDS 200).

For each sample 7 PCR tubes are analyzed. Each tube contains primers for the amplification of the internal control and for the target allele with the corresponding specific probes.

For each sample analyzed, individual alleles are detected by the kit and the corresponding serologic haplotype.

Celiac Gene Alleles identifies the presence of the alleles that codify for the DQ2 and DQ8 heterodimers that determine the risk of developing Celiac Disease. In addition the kit identifies the DQB1*02 homozygous status.

The **DQ2 heterodimer** is identified by the presence of the DQA1*05 and DQB1*02 allele. These alleles are commonly found in the DRB1*03-DQA1*05-DQB1*02 haplotype (cis-combination) or in the DRB1*11/12-DQA1*05-DQB1*0301/DRB1*07-DQA1*0201-DQB1*02 haplotypes (trans-combination).

The **DQ8 heterodimer** is identified by the presence of the DQA1*03 and DQB1*0302/0305 allele. These alleles are found in the DRB1*04-DQA1*03-DQB1*0302/0305 haplotype (cis-combination).

The presence of the DQB1*02 allele with a DQA1 allele diverse than DQA1*05 allele, identifies a DQ2 heterodimer different from primary associations with Celiac Disease. This DQ2 heterodimer, which determines a moderate risk of developing the disease, is found in the DRB1*07-DQA1*0201-DQB1*02 haplotype.

The risk is calculated by considering only the DQ primary associations with the disease.

VALIDATION

The test is not valid when the phrase "REPEAT ANALYSIS" appears on the BioRun touchscreen and the sample must be retested starting from the DNA extraction.

If you run the "PROCESS CONTROL" the word "POSITIVE" must appear on the display. If the word "NEGATIVE" appears, the test must be repeated for all the samples starting from the DNA extraction.

If you run the "CONTAMINATION CONTROL", results must indicate "NO CONTAMINATION".

If "CONTAMINATION" appears, the work session is considered not valid and the work area must be sanitized before re-beginning the procedure.

TROUBLESHOOTING

PROBLEM	CAUSE	SOLUTION
The words "REPEAT ANALYSIS" appear on the touchscreen	DNA absent	Repeat DNA extraction. If the problem continues repeat blood withdrawal
	Blood collected in the wrong tube (i.e. collected in a tube containing heparin)	Repeat blood sample withdrawal and store the sample in an EDTA-tube
	Agglutinated or coagulated blood sample	Do not use; repeat blood sample withdrawal
	Instruments (thermal cycler) not calibrated	Use and/or calibrate the lab instruments used for performing the tests
	PCR session failed due to improper use of the thermal cycler	Place PCR tubes directly onto thermal cycler plate, do not use PCR tube rack. Control amplification program
	Evaporation in reaction tubes, due to water loss and change in the concentration of the PCR reagents	Tighten and secure tubes with caps - vortex and centrifuge to collect drops before using thermal cycler
An error code appears on the touchscreen	The fluorescence values for the identification of the internal control and for the identification of the target allele are low. The primer mix located on the bottom of the tube has not been sufficiently dissolved	Repeat the test and vortex each tube of the strip adequately
	The fluorescence value for the identification of the target allele is in the uncertain range. The primer mix located on the bottom of the tube has not been sufficiently dissolved	Repeat the test and vortex each tube of the strip adequately
The word "NEGATIVE" appears in the process control	PCR session failed	You must repeat the test, starting from DNA lysis
The word "CONTAMINATION" appears in the contamination control	Contamination due to amplified products from previous work sessions	Sanitize the work area and instruments

Note: You can contact BioDiagene support team to receive more information.

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Baggrundsoplysninger:

EliA journal: Mummert, E., EliA Celikey – the first fully automated tTg-antibody test.

Andersen, R. et al. Screening for Coeliac Disease: Integration of technology and stakeholders.

Collin, P. 2004. Antibodies and Biopsy in the Diagnosis of Coeliac Disease. Lessons from the Celikey European Multicentre Study.

Wong, R.C.W. et al., 2002, A comparison of 13 guinea pig and human anti-tissue transglutaminase antibody ELISA kits, J. Clin. Pathol., 55, 488-494.

Williams L. et al., 2004. Guideline for the Diagnosis and Treatment of Celiac Disease in Children: Recommendations of the North American Society for Pediatric Gastroenterology, Hepatology and Nutrition.

Klinisk Info, 2006. Cøliaki (non tropisk sprue) (Klinisk Info. (tidl. Århus amt).

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Navigation buttons: Left arrow, Up arrow, Down arrow, Right arrow, Back, Information, Message.

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The image shows a software interface titled "TenKey". It features a "Test:" label above a text input field containing "Pipetting Depth". Below this, there are two input fields: "Lower limit:" with the value "0.00" and "Upper limit:" with the value "105.00". A large, empty input field is positioned below these, with the number "1.0" visible on its right side. At the bottom of the window is a numeric keypad with buttons for digits 0-9, a decimal point, and left and right arrow keys. To the right of the keypad are two buttons labeled "Cancel" and "Confirm".

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3. Vælg **Pipette Dispensing Test**.

Compensation value		Compensation value			
Sample	10123	Stop Solution			
Diluent	10000	Use			
Dilution ratio		<input checked="" type="radio"/> Sample	<input type="radio"/> Blank	<input type="radio"/> FluoroC	
<input checked="" type="radio"/> 1(Pre-Dil)	<input type="radio"/> 5	1:1	9963	1:3-1st	9913
<input type="radio"/> 2	<input type="radio"/> 10	1:2-1st	9912	1:3-2nd	10145
<input type="checkbox"/> Diluted Sample		1:2-2nd	10195	1:3-3rd	10145
Strip	10032	Reaction mixture (ELISA)	0		
Conjugate	10070	Development	10130		
Replicate	5				

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
Meddelelse 2-11

Message							
Previous	Next	Pipetting Test	Clogging Test		Stop		Mute
Select test (2-11)							
<p>Select "Pipetting test" or "clogging test".</p> <p>Pipetting Test: The Pipetting test is started.</p> <p>Clogging Test: The Clogging test is started.</p> <p>Stop: The following processing is discontinued and screen is returned to the menu.</p> <p>Mute: The warning sound is stopped.</p>							<p>Time</p> <p>10:21:54</p>
							<p>PrevPage</p> <p>NextPage</p>
					Information	Message	

Vælg Pipetting Test.

Meddelelse 1-14

Message							
Previous	Next		Start		Stop		Mute
Please set ImmunoCAP. (1-14) Number of CAP: 5							
<p>Please put it from the before one hole of conjugate wash1 position when exhaling on the immuno reaction wheel.</p> <p>Please put it from the CAP transfer position when exhaling on the enzyme reaction wheel.</p> <p>Two or more selections: Please put Cap sequentially.</p> <p>Immuno</p> <ol style="list-style-type: none"> 1. Sample 2. Diluent 3. Diluted Sample 4. Strip 5. Conjugate <p>Enzyme</p> <ol style="list-style-type: none"> 1. Development 2. Stop Solution or Reaction mixture in EIA <p>(note) Please use Cap where Development was pipetted when you do Reaction mixture in EIA. And, please measure the amount of the decrease of Cap before and after the examination.</p> <p>Start: The test is started.</p> <p>Stop: The following processing is discontinued and screen is returned to the menu.</p>							<p>Time</p> <p>10:22:30</p>
							<p>PrevPage</p> <p>NextPage</p>
					Information	Message	

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$$\text{Total CV (\%)}_{\text{within runs}} = 100 \frac{\sqrt{V_{Within}}}{M_{Tot}}$$

$$\text{Total (\%)} = 100 \frac{\sqrt{V_{Mean} + 0.5 \times V_{Within}}}{M_{Tot}}$$

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$$\text{Total CV (\%)}_{\text{within runs}} = 100 \frac{\sqrt{V_{\text{Within}}}}{M_{\text{Tot}}}$$

$$\text{Total (\%)} = 100 \frac{\sqrt{V_{\text{Mean}} + 0.75 \times V_{\text{Within}}}}{M_{\text{Tot}}}$$

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



Testprincip

 <p>The diagram shows a curved surface labeled 'ImmunoCAP' with red dots representing allergens. Blue Y-shaped molecules labeled 'IgE Antibody' are shown binding to these allergens.</p>	<p>Det allergen der undersøges for, og som er kovalent koblet til celluloserøret i ImmunoCAP, reagerer med det specifikke IgE i patientprøven.</p>
 <p>The diagram shows the 'ImmunoCAP' support with allergens. Some blue Y-shaped molecules are bound to the allergens, while others are free in the solution, representing non-specific IgE being washed away.</p>	<p>Efter at ikke-specifikt IgE er vasket væk, tilsættes enzymmarkerede antistoffer mod IgE for at danne et kompleks.</p>
 <p>The diagram shows the 'ImmunoCAP' support with allergens. Blue Y-shaped molecules are bound to the allergens, and purple Y-shaped molecules (enzyme-labeled anti-IgE) are shown binding to the blue molecules.</p>	<p>Efter inkubationen vaskes ubundet enzym-anti-IgE væk, og det bundne kompleks inkuberes derefter med et udviklingsmiddel.</p>
 <p>The diagram shows the 'ImmunoCAP' support with allergens. Blue Y-shaped molecules are bound to the allergens, and purple Y-shaped molecules are bound to the blue molecules. Yellow starburst shapes are shown near the bound complex, indicating a visible reaction.</p>	<p>Når reaktionen er stoppet, måles eluatets fluorescens. Jo højere fluorescens, desto mere specifikt IgE er der i prøven. For evaluering af testresultaterne sammenlignes responsen for patientprøverne direkte med responsen for kalibratorerne.</p>

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
1. Fejlliste
2. Systemparameterindstillinger
3. Temperaturovervågning

Rediger doser

The image shows a software dialog box titled "Carrier Barcode". It has a light blue background and a dark blue title bar with a close button (X) in the top right corner. Inside the dialog, there are two text input fields. The first is labeled "ImmunoCAP Carrier ID" and contains the text "8819901". The second is labeled "Amount:" and contains the text "16". Below these fields is a radio button that is currently selected, with the label "Barcode Reading" to its right. At the bottom of the dialog, there are two buttons: "OK" on the left and "Cancel" on the right.

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
Prøveoplysninger

The screenshot shows a software window titled "Sample Information". Inside, there are two input fields: "Sample Rack Position:" and "Sample Rack ID:". The "Sample Rack ID:" field contains the value "0003". Below these fields is a table with the following data:

No.	ID	Test name	Status	Dilution factor	Priority
1	R7762008	t2	Completed		Normal
2	R7762008	t3	Completed		Normal
3	R7762008	a-IgE	Completed		Normal
4	R7762009	t1	Completed		Normal
5	R7762009	t2	Completed		Normal
6	R7762009	t3	Completed		Normal
7	R7762009	a-IgE	Completed		Normal
8	R7762010	t2	Completed		Normal
9	R7762010	a-IgE	Completed		Normal

At the bottom of the window, there are three buttons: "Back", "Information", and "Message".

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Oplysninger om opbevaring af ImmunoCAP/EliA Well-rør

ImmunoCAP Carrier Storage Tray info.

Tray Position: Tray ID:

Pos	Status	Testname	Amount	Lot number	Lot specific	Expiration Id.	Use	Days loaded
1	Empty	a_IgE	0	81801		2007/11/01		8
2	OK	I1	16	36101		2007/11/01		8
3	OK	a_IgE	16	86501		2007/11/01		8
4	OK	a_IgE	16	86501		2007/11/01		8
5	OK	Gcal	12	BFN01		2006/11/28		8
6	OK	a-IgE	7	86301		2007/11/01		8
7	OK	I3	12	40401		2007/11/01		8
8	OK	w1	9	40701		2007/11/01		8
9	OK	a_IgE	16	86501		2007/11/01		8
10	OK	I1	5	406ZZ		2088/01/01		8
11	OK	I2	13	392ZZ		2088/01/01		8
12	OK	I4	16	393ZZ		2088/01/01		8

Navigation buttons: Left arrow, Up arrow, Right arrow, Down arrow, Magnifying glass.

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Fejlliste

Error List			
Date Time	Error/Warn	Error No.	Information
2006/11/07 13:00:55	Error	00001191	The difference of volume : (0-1191) : 4 Position:2
2006/11/07 12:58:23	Error	00002187	Barcode is abnormal - Barcode is not recognized by this screen or instrumen
2006/10/31 13:23:47	Error	00002052	A rack in process was wrongly removed from the sample load area.(0-2052) I
2006/10/31 12:31:16	Error	00002052	A rack in process was wrongly removed from the sample load area.(0-2052) I
2006/10/31 12:30:31	Error	00001191	The difference of volume : (0-1191) : 9 Position:1
2006/10/31 11:06:13	Error	00002052	A rack in process was wrongly removed from the sample load area.(0-2052) I
2006/10/30 15:35:03	Error	00002187	Cod.Barras anormal - Cod.Barras no es reconocido por el instrumento. (0-2187
2006/10/30 15:11:21	Error	00002060	A rack with the same ID is already loaded. This rack is passed.(0-2060) Rack
2006/10/30 15:10:39	Error	00002060	A rack with the same ID is already loaded. This rack is passed.(0-2060) Rack
2006/10/30 14:37:50	Error	00002187	Barcode is abnormal - Barcode is not recognized by this screen or instrumen
2006/10/30 12:59:34	Error	00002187	Barcode is abnormal - Barcode is not recognized by this screen or instrumen
2006/10/30 12:50:22	Error	00002178	There is no free position of a strip bottle. Please take out a strip bottle (0-2178
2006/10/30 12:39:03	Error	00002052	A rack in process was wrongly removed from the sample load area.(0-2052) I
2006/10/30 12:38:56	Error	00002052	A rack in process was wrongly removed from the sample load area.(0-2052) I
2006/10/30 12:38:48	Error	00002052	A rack in process was wrongly removed from the sample load area.(0-2052) I
2006/10/30 12:38:37	Error	00002052	A rack in process was wrongly removed from the sample load area.(0-2052) I

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Brugervedligeholdelse

The screenshot shows the 'User Maintenance' interface with the following sections:

- Prime:** An illustration of a blue cylindrical component with two purple clamps.
- Blank Measurement:** A red-bordered gauge with a needle pointing to approximately 10 on a scale from 0 to 20.
- Daily Rinse:** A calendar grid for the month of June. The 1st, 2nd, 3rd, 4th, 10th, 11th, 12th, 13th, 14th, 15th, 16th, 17th, 18th, 19th, 20th, 21st, 22nd, 23rd, 24th, 25th, 26th, 27th, 28th, 29th, 30th, and 31st are highlighted in blue.
- Weekly Rinse:** A calendar grid for the month of June. The 1st, 2nd, 3rd, 4th, 10th, 11th, 12th, 13th, 14th, 15th, 16th, 17th, 18th, 19th, 20th, 21st, 22nd, 23rd, 24th, 25th, 26th, 27th, 28th, 29th, 30th, and 31st are highlighted in blue.
- Monthly Maintenance:** A calendar grid for the month of June. The 1st, 2nd, 3rd, 4th, 10th, 11th, 12th, 13th, 14th, 15th, 16th, 17th, 18th, 19th, 20th, 21st, 22nd, 23rd, 24th, 25th, 26th, 27th, 28th, 29th, 30th, and 31st are highlighted in blue.

At the bottom of the screen, there are three buttons: 'Back', 'Information', and 'Message'.

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- PRIME
- BLANK MEASUREMENT
- DAILY RINSE
- WEEKLY RINSE
- MONTHLY MAINTENANCE

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FluoroC-kørsel

Status	Lot number	Expiration Id...	State	
In use	95201	2007/11/01		

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FluoroC-måling

The screenshot shows the 'FluoroC Measurement' window. It contains several input fields and buttons for configuring and executing the measurement. At the top, there are fields for 'Rinse Blank' (0.13) and 'Upper' (1.00), each with an 'OK' button. Below these is a section titled 'FluoroC values' which contains fields for 'Replicate(1)' (473.47), 'Replicate(2)' (469.51), 'Replicate(3)' (470.42), 'CV%' (0.44), 'RepDeviation' (-5.58), 'InnerLimit' (4.00), 'OuterLimit' (10.00), and 'Target' (D6). Each of these fields has an associated 'OK' button. At the bottom of the window is a row of buttons: 'Pause', 'Abort', a small icon button, 'Back', 'Information', and 'Message'.

Field	Value	Action
Rinse Blank	0.13	OK
Upper	1.00	
Replicate(1)	473.47	
Replicate(2)	469.51	
Replicate(3)	470.42	
CV%	0.44	OK
RepDeviation	-5.58	
InnerLimit	4.00	NG
OuterLimit	10.00	OK
Target	D6	

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
Valg af datoformat

The screenshot shows a 'Select Data' dialog box with a title bar. Inside, there is a 'Data Name:' label followed by a text box containing 'Date Format'. Below this is a 'Select Item' section with a light blue background. It contains three radio button options for date formats: 'yyyy/mm/dd(yy/mm/dd)' (selected), 'dd/mm/yyyy(dd/mm/yy)', and 'mm/dd/yyyy(mm/dd/yy)'. At the bottom of the dialog is a bar with four buttons: 'Back', 'Information', 'Message', and a small icon button.

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
Måleenhed for temperatur

The image shows a software window titled "Select Data". Inside, there is a text field labeled "Data Name:" containing the text "Temperature Unit". Below this is a larger area labeled "Select Item" which contains two radio button options: "Celsius" (which is selected, indicated by a filled circle) and "Fahrenheit" (which is unselected, indicated by an empty circle). At the bottom of the window, there is a row of three buttons: "Back", "Information", and "Message". To the left of these buttons are three empty rectangular slots and a small icon of a triangle pointing up.

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Rørindstillinger

The screenshot shows a software window titled "Parameter - Tube Setting". It contains two columns of input fields: "Bottom Level[0.0-105.0mm]" and "Pipetting Depth[mm]". The rows are labeled as follows:

	Bottom Level[0.0-105.0mm]	Pipetting Depth[mm]
Tube 1 (normal)	1.1	1.0
Tube 2 (pediatric)	39.0	1.5
Tube 3	7.8	1.1
Tube 4	1.1	1.0
Tube 5 (reserved)	1.1	1.0
QC vial 1	7.8	1.1
QC vial 2 (small)	39.0	104.0

At the bottom of the window, there are four buttons: "Register", a button with a triangle icon, "Back", "Information", and "Message".

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
The screenshot shows a software window titled "Select Multiple Item". Inside, there is a text field labeled "Data Name:" containing the text "After Assay: Washing/Rinse/Soak". Below this is a larger box labeled "Select Item" which contains three items, each with a checked checkbox: "Washing", "Rinse", and "Soak". At the bottom of the window, there are three buttons: "Back", "Information", and "MESSAGE".

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
The screenshot shows a 'Select Data' dialog box. At the top, there is a 'Data Name:' label followed by a text box containing 'Accept low frequent test from IDM'. Below this is a 'Select Item' label followed by a large rectangular area. Inside this area, there are two radio buttons: 'No' on the left and 'Yes' on the right. The 'Yes' radio button is selected. At the bottom of the dialog box, there is a row of buttons: 'Back', 'Information', and 'MESSAGE'. To the left of these buttons are three empty rectangular boxes and a small icon of a computer monitor.

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
Select Data

Data Name: After Assay: Automatic Power Off Instrument

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The screenshot shows a software window titled "Select Data". Inside, there is a text field labeled "Data Name:" containing the text "Eject Empty Carriers to Waste Box". Below this is a section labeled "Select Item" which contains two radio buttons: "No" (which is selected) and "Yes". At the bottom of the window, there is a row of three buttons: "Back", "Information", and "Message".

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Time set less than 30 minutes from actual time refers to start up next week!

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Lagersektionsliste for ImmunoCAP/EliA Well

Carrier Storage Tray List	
Carrier Storage Tray List	
Tray ID	ImmunoCAP Group No
TA	0,1,2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19
TB	0,1,2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19
TC	0,1,2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19
TD	0,1,2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19
TE	0,1,2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19
TF	0,1,2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19
TG	0,1,2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19
TH	0,1,2,3,4,5,6,7,8,9,10,11,12,13,14,15,16,17,18,19
TI	20,21,22,23,24,25,26,27,28,29,30,31,32,33,34,35,36,37,38,39
TJ	20,21,22,23,24,25,26,27,28,29,30,31,32,33,34,35,36,37,38,39

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
Error Modify

Err No.

Comment
There is no response from the subsystem. (0-1001) Item:%d

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The image shows a 'Select Data' dialog box. At the top, there is a label 'Data Name:' followed by a text input field containing 'UPS use'. Below this is a section titled 'Select Item' which contains two radio buttons: 'No' (which is selected) and 'Use'. At the bottom of the dialog, there are three buttons: 'Back', 'Information', and 'Message'.

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
- Heater On/Off (Standby Mode)
- Heater Off Timer

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The screenshot shows the 'Assay Processing' window. At the top, 'Status' is set to 'Running' and 'Active Technology' is 'ImmunoCAP'. Below is a table with 10 rows of assay steps. At the bottom, there are navigation arrows and buttons for 'Load', 'End assay', 'Change Technology', 'Back', 'Information', and 'Message'.

No.	ID	Method	Test	Dilution	Time	Process Status
1	CAL-0.35	slgE	a_lgE		12:55	Sample pipetting
2	CAL-0.35	slgE	a_lgE		12:55	Sample pipetting
3	CAL-0.70	slgE	a_lgE		12:55	Sample pipetting
4	CAL-0.70	slgE	a_lgE		12:55	Pre wash
5	CAL-3.50	slgE	a_lgE		12:55	Pre wash
6	CAL-3.50	slgE	a_lgE		12:55	Pre wash
7	CAL-17.5	slgE	a_lgE		12:55	ImmunoCAP/EliA Well dispense
8	CAL-17.5	slgE	a_lgE		12:55	ImmunoCAP/EliA Well dispense
9	CAL-50.0	slgE	a_lgE		12:55	ImmunoCAP/EliA Well dispense
10	CAL-50.0	slgE	a_lgE		12:55	ImmunoCAP/EliA Well dispense

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The screenshot shows the 'ASSAY PROCESSING PAUSE/CONTINUE' dialog box. It has two sections: 'End assay choice' and 'Rinse option'. The 'End assay choice' section has three buttons: 'End after all loaded samples are processed' (highlighted), 'End after dispensed ImmunoCAP/EliA Well are processed', and 'Abort'. The 'Rinse option' section has two buttons: 'Daily rinse' (highlighted) and 'Primed (Standby)'. At the bottom are 'OK' and 'Cancel' buttons.

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Quality Select Data IC250 (ImmunoCAP 250) E-I-G Print Back ESC

Calibrators **Curve Controls** **Quality Controls** **Quality Club**

Period/Sample No	Method name
0806	sgE
0807	E-I-G
0807	Tryp
0808	ECF
0808	sgE
200809	E-I-G

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Request List Find Print Back ESC

Samples **Requestors**

Sample ID	Status	Sample Date	Type	Tube	Pre-dil	Rack ID	Position	Error Occured	Request ID
00001	Loaded	2008-06-12	Patient	Normal	1			No	00001
00001	Not started	2008-05-10	Patient	Normal	1			No	00001
00002	Not started	2008-05-28	Patient	Normal	1			No	00002
0806-225656	Not started	2008-08-29	Quality club	Normal	1			No	
0806-258562	Not started	2008-08-29	Quality club	Normal	1			No	
0807-245624	Not started	2008-08-29	Quality club	Normal	1			No	
0807-255686	Not started	2008-08-29	Quality club	Normal	1			No	
0807-256354	Not started	2008-08-29	Quality club	Normal	1			No	
0808-45365	Not started	2008-08-29	Quality club	Normal	1			No	
0808-532677	Not started	2008-08-29	Quality club	Normal	1			No	
0808-532678	Approved	2008-08-29	Quality club	Normal	1	0001	1	No	
121324	Not started	2008-09-02	Patient	Normal	1			No	121324
123123	Not started	2008-09-02	Patient	Normal	1	0001	2	No	123123
123124	Not started	2008-09-02	Patient	Normal	1	0001	3	No	123124
200809-1234	Not started	2008-09-09	Quality club	Normal	1			No	

New **Open** **Delete**

☒ All
☐ Not started
☐ Processing
☐ Ready
☐ Approved
☐ Rejected
☐ Reported
☐ Reflex

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Date

Origin

☐ High priority
☐ Hospitalized
☐ Available in laboratory

Save

Tests

Result graph

Test	Method	Status	Conc	Class	Quotient	Cut-off	Cut-off 2	Repl	Dil	Full Name	Reflex
sy	EI-G	Not started	—	—				1	100	EIA Symphony	No
ro	EI-G	Not started	—	—				1	100	EIA Ro	No
la	EI-G	Not started	—	—				1	100	EIA La	No

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Print

ESC Back

	Sample Id	Pediatric	Request	Skip barcode reading
1	200809-123456	<input type="checkbox"/>		<input checked="" type="checkbox"/>
2		<input type="checkbox"/>		<input type="checkbox"/>
3		<input type="checkbox"/>		<input type="checkbox"/>
4		<input type="checkbox"/>		<input type="checkbox"/>
5		<input type="checkbox"/>		<input type="checkbox"/>
6		<input type="checkbox"/>		<input type="checkbox"/>
7		<input type="checkbox"/>		<input type="checkbox"/>
8		<input type="checkbox"/>		<input type="checkbox"/>
9		<input type="checkbox"/>		<input type="checkbox"/>
10		<input type="checkbox"/>		<input type="checkbox"/>

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The screenshot shows the 'Quality Club' software interface. At the top, there is a 'Select Data' dropdown menu set to 'UC1000#1 (ImmunoCAP 1C)' and a 'Select method' dropdown menu. To the right are 'Print' and 'Back' buttons. Below these are four tabs: 'Calibrators', 'Curve Controls', 'Quality Controls', and 'Quality Club'. The 'Quality Club' tab is selected. It contains a table with two columns: 'Period/Sample No' and 'Method name'. The table lists five entries: 200801-1, 200802-1, 200803-1, 200804-1, and 200808, all with method names 'slgE' or 'EHG'. To the right of this table are 'New' and 'Delete' buttons. Below the table is another table with columns: 'SampleId', 'Sample Lot No', 'Test', 'Result', 'Quotient', 'Instrument', 'InstrumentType', and 'SerialNumber'. This table lists three entries for sample '200808-BJ386' with test results 'ce', 'dn', and 'la', all showing a result of '4 U/ml'. To the right of this table is a 'Details' button.

Period/Sample No	Method name
200801-1	slgE
200802-1	slgE
200803-1	slgE
200804-1	slgE
200808	EHG

SampleId	Sample Lot No	Test	Result	Quotient	Instrument	InstrumentType	SerialNumber
200808-BJ386	BJ386	ce	4 U/ml		ImmunoCAP 1000	UC1000#1	
200808-BJ386	BJ386	dn	4 U/ml		ImmunoCAP 1000	UC1000#1	
200808-BJ386	BJ386	la	4 U/ml		ImmunoCAP 1000	UC1000#1	

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Other reagents used

Calibrator Kit/Strip Lot No.	A558A	A558A	A558A
Curve Control Kit/Strip Lot No.	APVCA	APVCA	APVCA
Calibrator ImmunoCAP/Well Lot No.	R6522	R6522	R6522
Conjugate Lot No.	R0832	R0832	R0832

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311

Quality Club			
ELIA IgG			
Laboratory code	200808		
Period/Sample No.	UC1000*1		
Instrument	ImmunoCAP 1000		
Type	ImmunoCAP 1000		
SerialNumber			
Results			
Sample Lot No.	BJ386	BJ386	BJ386
Test (ImmunoCAP/Allergen/Antigen)	dn	ce	la
ImmunoCAP/Elia Well Lot No.			
Results	4 U/ml	4 U/ml	4 U/ml
Other reagents used			
Calibrator Kit/Strip Lot No.	A5SBB	A5SBB	
Curve Control Kit/Strip Lot No.	A0VAA	A0VAA	
Calibrator ImmunoCAP/Well Lot No.	B65BS	B65BS	
Conjugate Lot No.	BGRCP	BGRCP	

Fax	Fax number not defined		
			

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
Identity	Test	Measured	Approval	Dilution	Resp	Conc	Class	Cut-off	Quot	Cut-Off2	Error
gbm h	gb	16:35:05	Not appro...	100	2630	53 U/ml	Positive				OK
DNA L	dn	16:48:04	Not appro...	10	316	3.3 U/ml	Negati...				OK
DNA	dn	16:49:03	Not appro...	10	2230	26 U/ml	Positive				OK
DNA	dn	16:50:03	Not appro...	10	4326	53 U/ml	Positive				OK
CCPL	cp	16:51:03	Not appro...	100	530	4.0 U/ml	Negati...				OK
CCPL	cp	16:52:03	Not appro...	100	2237	18 U/ml	Positive				OK
CCPL	cp	16:53:03	Not appro...	100	6708	61 U/ml	Positive				OK
GGL L	Ggl	17:08:01	Not appro...	100	845	4.9 U/ml	Negati...				OK
GGL	Ggl	17:09:01	Not appro...	100	3946	25 U/ml	Positive				OK
GGL	Ggl	17:10:01	Not appro...	100	5698	38 U/ml	Positive				OK
MPO	mp	17:11:01	Not appro...	50	206	3.3 U/ml	Negati...				OK
MPO	mp	17:12:01	Not appro...	50	2006	36 U/ml	Positive				OK
MPO	mp	17:13:01	Not appro...	50	3455	65 U/ml	Positive				OK
PR3 L	pr	17:14:01	Not appro...	50	353	5.2 U/ml	Negati...				OK
PR3	pr	17:15:01	Not appro...	50	2568	42 U/ml	Positive				OK
GBM	gb	17:16:00	Not appro...	100	365	6.8 U/ml	Negati...				OK
GBM	gb	17:22:00	Not appro...	100	1676	33 U/ml	Positive				OK
GCY	Gcy	17:23:00	Not appro...	100	2908	58 U/ml	Positive				OK
GBM	gb	17:27:59	Not appro...	100	2535	51 U/ml	Positive				OK
GBM	gb	17:28:59	Not appro...	100	2701	55 U/ml	Positive				OK
GBM	gb	17:29:59	Not appro...	100	2526	51 U/ml	Positive				OK
GBM	gb	17:30:58	Not appro...	100	2679	54 U/ml	Positive				OK
GBM	gb	17:31:58	Not appro...	100	2710	55 U/ml	Positive				OK
GBM	gb	17:32:58	Not appro...	100	2731	55 U/ml	Positive				OK
GBM	gb	17:33:58	Not appro...	100	2566	52 U/ml	Positive				OK
GBM	gb	17:34:58	Not appro...	100	2681	54 U/ml	Positive				OK
GBM	gb	17:35:58	Not appro...	100	2844	58 U/ml	Positive				OK
PR3 S	pr	17:45:57	Not appro...	50	3314	55 U/ml	Positive				OK
PR3 S	pr	17:46:57	Not appro...	50	1662	26 U/ml	Positive				OK
PR3 S	pr	17:47:57	Not appro...	50	1172	18 U/ml	Positive				OK

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Phadia AB er indehaver af ophavsretten til nærværende manual, © Phadia AB, alle rettigheder forbeholdes. Intet materiale må reproducere eller genudgives uden skriftlig tilladelse fra Phadia AB.